



eFOLDi<sup>®</sup>

## eFOLDi Explorer

eFOLDi Explorer is one of the lightest, most compact, folding mobility scooter in its class, weighting only 17kg (without battery) and no larger than a small suitcase when folded. The eFOLDi Explorer is fully UK road legal with a top speed of 8mph and a range of 14 miles. eFOLDi Explorer fits easily into the smallest car and is fully airline compliant - the perfect travel companion.

**Enjoy freedom, friendship and fun!**



eFOLDi<sup>®</sup>

### Features

Speed	4mph (6kmh) / 8mph (12kmh)
Range	14 miles / 22 Km
Carry capacity	120Kg / 265 lbs / 19 Stone
Climbing ability	1:6 Gradient (9.5 degrees)
Folded dimensions	66 x 56 x 40 cm / 26 x 22 x 15.6 inches
Tyre type	Pneumatic
Battery standard	12AH Li-ion
Motor	240W 24V Hub Motor
Ignition	Key operated